

PREMIUM WELLNESS FOOTWEAR

CLASS 1 MEDICAL DEVICE
PATENTED COOLMAX
CARBON LINING

DESIGNED FOR

- TYPE 1 & 2 Diabetes Mellitus
 - Rheumatoid
 - Bone Diseases
 - Digital Deformties
 - Hallux Abductu Valgus
 - Hammer Toes
 - Claw Toes
 - Predislocations
 - Bone coalitions

IN STOCK PROGRAM
BY THE PAIR





COOLMAX CARBON

We have developed an innovative footwear collection, which combines functional footwear designs for various types of pathologies with construction materials.

Our current collection uses COOLMAX® fabric lining to which we add Carbon threads.

COOLMAX® polyester fibers have long been known for their high breathability due to their hollow fiber design and aeration channels, which helps to release moisture quickly and efficiently.

By combining COOLMAX® fibers with Carbon fibers we have innovated a one of a kind fabric with the ability to create a dry, and airy environment while maintaining freshness.

In a test for Staphylococcus aureus and Klebsiella pneumoniae, two bacteria that can proliferate under normal conditions of humidity and temperature caused by sweating our Coolmax- Carbon fabric has demonstrated bacteriostatic properties and maintained its freshness.

Unlike other fabrics, Coolmax-Carbon fabric does not receive any chemical treatments such as microencapsulation, or ion application. These treatments degrade over time resulting in the fabric losing its properties.

Our innovative footwear collection that combines functional designs and advanced Coolmax-Carbon fabric ensures the general comfort and safety of the wearers.


HEAD OFFICE


MASCIOTRA DESIGN INC

T. 514-758-7463

INFO@MASCIOTRADESIGN.COM

DANIELE MASCIOTRA





Zephir

Description:

Zephir® footwear has been designed for persons with wide forefoot, who consequently need an easy way to put shoes on and a material that provides them comfort. It is also designed for big digital deformities.

Application:

The Zephir® footwear is a health use product recommended for persons with fitting problems in the forefoot. It is also designed with adaptable adjustments and special fabrics for systemic diseases, such as persons with Type I and II Diabetes Mellitus, rheumatoid or bone diseases and for its antibacterial properties and its special widths. It is especially indicated for digital deformities, such as; Hallux Abductus Valgus, hammer or claw toes, predislocations of metatarsophalanges and bone coalitions on the back of the foot.

Classification:

Zephir® footwear is a class I (rule 1) medical device. The evaluation of conformity was made following Annex VII (EC Declaration of Conformity), of the Council Directive 93/42/EEC and Directive 2007/47.





SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48
(a) FOOTWEAR LENGTH MM	240	247	253	260	267	273	280	287	293	300	307	313	320
(b) FOOTWEAR LENGTH INCHES	9,45	9,72	9,96	10,23	10,51	10,75	11,02	11,3	11,53	11,81	12,1	12,32	12,6



black







Description:

Ortho® footwear has been designed for persons with forefoot and toe problems who require a special width, an easy way to put them on and material that provides more comfort.

Application:

Ortho® footwear is a health use product recommended for persons with digital deformities, such as; hammer or claw toes, predislocated or dislocated toes, Hallux Abductus Valgus, Quintus Varus, deorsumducted or sursumducted toes and pronounced instep. It is also recommended for systemic diseases like Type I and II Diabetes Mellitus and rheumatoid or bone diseases, due to its antibacterial properties and the fit and adjustment provided by their Velcro and their special widths.

Classification:

Ortho® footwear is a class I (rule 1) medical device. The evaluation of conformity was made following Annex VII (EC Declaration of Conformity), of the Council Directive 93/42/EEC and Directive 2007/47.





SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48
(a) FOOTWEAR LENGTH MM	240	247	253	260	267	273	280	287	293	300	307	313	320
(b) FOOTWEAR LENGTH INCHES	9,45	9,72	9,96	10,23	10,51	10,75	11,02	11,3	11,53	11,81	12,1	12,32	12,6





Description:

Wallaby® footwear has been designed for persons with various walking functionality problems and who require a special width and material that provides them greater comfort.

Application:

Wallaby® footwear is a health use product for persons with Type I and Type II Diabetes Mellitus; for its anti-bacterial properties, special widths and technical linings. For these reasons, it is recommended for persons with a deformity in the forefoot or toes and systemic diseases that affect the foot, such as rheumatoid or bone diseases.

Classification:

Wallaby® footwear is a class I (rule)1 medical device. The evaluation of conformity was made following Annex VII (EC Declaration of Conformity), of the Council Directive 93/42/EEC and Directive 2007/47.





SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48
(a) FOOTWEAR LENGTH MM	240	247	253	260	267	273	280	287	293	300	307	313	320
(b) FOOTWEAR LENGTH INCHES	9,45	9,72	9,96	10,23	10,51	10,75	11,02	11,3	11,53	11,81	12,1	12,32	12,6



CH31

Description:

The CH31® footwear has been designed for persons with problems of ulcerations or ischemia, which require a special width and material that provides a better fit and comfort.

Application:

The CH31® footwear is a health use product for persons with Type I and Type II Diabetes mellitus, known for its antibacterial properties, special widths and easy

They are also suitable for persons with systemic diseases with feet deformities, such as; rheumatic bone and circulation disorders. Their Velcro in the heel area facilitates their fitting and it is also indicated for pathologies in this zone, such as; Haglund's disease, lesions on the retro Achilles bursa and the calcaneal bone. They may also be used for problems related to trauma accidents that require special widths and a very comfortable adjustment and fit.

Classification:

The CH31® footwear is a class I (rule 1) medical device. The evaluation of conformity was made following Annex VII (EC Declaration of Conformity), of the Council Directive 93/42/EEC and Directive 2007/47.





SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48
(a) FOOTWEAR LENGTH MM	240	247	253	260	267	273	280	287	293	300	307	313	320
(b) FOOTWEAR LENGTH INCHES	9,45	9,72	9,96	10,23	10,51	10,75	11,02	11,3	11,53	11,81	12,1	12,32	12,6

This style has been designed to adjust to fit widths from ww (2W) through wwwwww (6W).







AVAILABLE IN STOCK IN MONTREAL BY THE PAIR

Marc

Description:

the Marc® footwear has been designed for persons with problems of width in the forefoot and problems in the heel area. It is also designed for easy fitting and more comfort.

Application:

The Marc® footwear is a health use product recommended for persons with problems of fitting in the forefoot. It is also designed with adaptable adjustments and special fabrics for systemic diseases for persons with Type I and II Diabetes Mellitus, rheumatoid, bone or circulation diseases and for its antibacterial properties. It is especially indicated for enlargements on the heels, such as Haglung's disease and calcaneal calcifications.

Classification:

Marc® footwear is a class I (rule 1) medical device. The evaluation of conformity was made following Annex VII (EC Declaration of Conformity), of the Council Directive 93/42/EEC and Directive 2007/47.





SIZE	36	37	38	39	40	41	42	43	44	45	46	47	48
(a) FOOTWEAR LENGTH MM	240	247	253	260	267	273	280	287	293	300	307	313	320
(b) FOOTWEAR LENGTH INCHES	9,45	9,72	9,96	10,23	10,51	10,75	11,02	11,3	11,53	11,81	12,1	12,32	12,6



black

